

# Imaging Mueller Matrix Ellipsometry for the Characterization of Micro-Structured Anisotropic Thin-Film Samples


Matthias Duwe<sup>1</sup>, Sebastian Funke<sup>1</sup>, Christian Röling<sup>1</sup>, Peter H. Thiesen<sup>1</sup>, Aday J. Molina-Mendoza<sup>2</sup>, Andres Castellanos-Gomez<sup>3</sup>, Anna Sinterhauf<sup>4</sup>, Martin Wenderoth<sup>4</sup>  
<sup>1</sup>Accurion GmbH, Stresemannstr. 30, 37079 Goettingen, Germany; <sup>2</sup>Universidad Autonoma de Madrid, Departamento de Fisica de la Materia Condensada, Campus Universitario de Cantoblanco, 28049 Madrid, Spain; <sup>3</sup>IMDEA Nanoscience, C/ Faraday 9, Campus Universitario de Cantoblanco, 28049 Madrid, Spain; <sup>4</sup>IV. Physical Institute, Georg-August-University of Goettingen, Friedrich-Hund-Platz 1, 37077 Goettingen, Germany

## Introduction

Thin-film layers of black phosphorus - the most stable allotrope of the element phosphorus - is a promising material for future semiconductor electronics based on novel 2D-materials. Ultra-thin films or even monolayers of black phosphorus (b-P) may be fabricated by mechanical exfoliation from bulk material (similar to graphene). 2D-b-P is of particular interest for semiconductor electronics as it features a widely tunable bandgap via the layer thickness (from ~2 eV for monolayers down to 0.3 eV for bulk material). Thus, it closes the energy gap between the other known 2D-materials graphene (zero bandgap) and transition-metal dichalcogenides (bandgap typically 1-2 eV) such as MoS<sub>2</sub>. [1]

Regarding its optical properties, b-P stands out from other 2D-materials as well because its crystal structure causes an in-plane optical anisotropy in the visible range of the electromagnetic spectrum. However, the ellipsometric characterization of exfoliated flakes of b-P is extremely challenging as flake sizes are typically on the order of some 10µm, and the thickness often varies within a single flake.

In this work, for the first time we report the ellipsometric characterization of micro-structured b-P-flakes by means of spectroscopic Imaging Mueller Matrix Ellipsometry (IMME). It features much higher x-y-resolution compared to conventional ellipsometers since for IMME this resolution is not limited by the size of the illuminating spot. By using IMME with a microscope objective lens, the measurement of anisotropic refractive indices and the orientation of the optical axes becomes feasible even on microstructured thin-film samples such as b-P-flakes.



Mueller-matrix map at  $\lambda=501$  nm (left),  $\lambda$ -spectra obtained from spectroscopic mapping and evaluation of highlighted regions (top). Sample: b-P multilayer flake.


## Sample description

**Black phosphorus:**  
 ✓ black phosphorus consists of sheets of phosphorus atoms forming a puckered honeycomb lattice  
 ✓ Interlayer-bonds are weak (Van-der-Waals bonding)  
 ✓ mechanical exfoliation produces thin-film flakes or even monolayers


**Sample fabrication:**  
 ✓ thin-film flakes of black phosphorus were prepared by mechanical exfoliation (peeling off) from bulk material applying a polydimethylsiloxane (PDMS) stamp [2]  
 ✓ transfer of flakes onto a SiO<sub>2</sub>/Si-substrate by pressing PDMS stamp onto the substrate followed by a slow peel-off

The samples were prepared by Aday J. Molina-Mendoza and Dr. Andres Castellanos-Gomez.

**Ellipsometric modeling:**  
 ✓ describe b-P as orthorhombic, biaxial layer  
 ✓ Parameter:  
 ✓ complex refractive indices for in-plane and out of plane components of index ellipsoid  
 ✓ orientation of in-plane principal axes (=offset of rotational positioning stage)  
 ✓ Remark: b-P layer thickness will be determined independently (c. AFM-measurements)




SiO<sub>2</sub> layer thickness 292 nm (obtained from spectroscopic IMME measurements)



Orientation of in-plane optical axes relative to s-polarization in ellipsometric measurement

## Imaging MM Ellipsometry

**Imaging Ellipsometry**  
 ✓ combination of ellipsometry and microscopy  
 ✓ illuminating side: monochromatic illumination (collimated or slightly focused, no microspot required)  
 ✓ detecting side: microscope objective lens in detection beam path images sample onto camera (CCD, CMOS)



For an introduction to our technology please visit our technology page



**Used setup**  
 ✓ imaging ellipsometer in PCSA-configuration (nanofilm\_ep4, Accurion)  
 ✓ objective lens: 50x magnification  
 ✓ automated x-y-θ-sample stage  
 ✓ light source: super-continuum laser (SuperK extreme EXR-20, NKT photonics) with accusto-optic tunable filter (SuperK select, NKT photonics,  $\lambda=460-700$  nm)





For additional information, please visit our website




**Imaging Mueller Matrix Ellipsometry (IMME):**  
 ✓ operation in rotating compensator mode  
 ✓ measurement of 3x4-Mueller Matrix images (maps of the first 3 rows of normalized MM, normalization  $m_{11}=1$ )

## Imaging Mueller Matrix Ellipsometry $\theta$ -scans

✓ non-zero off-diagonal blocks of MM-micrographs reveal the b-P-flake's anisotropy (wavelength of illumination  $\lambda=550$  nm)  
 ✓  $\theta$ -scan: rotate the flake around its center point & record MM-micrographs at each position (rotation from 0° to 340° in 20°-steps)  
 ✓ evaluating central region of interest (ROI, c. fig. right) on each MM-micrograph for all  $\theta$ -positions yields Mueller-Matrix  $\theta$ -spectrum  
 ✓ observed periodicity of 180° complies with an orthorhombic, biaxial model for b-P (one axis orthogonal to the sample surface)  
 ✓ using the flake thickness obtained from AFM-measurements (43 nm) for optical modelling yields the in-plane refractive indices and the orientation of the in-plane principal axes of the index ellipsoid



Micromaps of normalized MM ( $m_{11}=1$ ) at  $\theta=0^\circ$  and  $\lambda=550$  nm. Mind different scale bars for block diagonal (blue frames) and block off-diagonal (green frames) matrix elements. Arrows in  $m_{34}$ -graph indicate the orientation of in-plane optical axes obtained from ellipsometric modelling of  $\theta$ -scan.



$\theta$ -spectrum of the normalized 3x4 MM of black phosphorus at different angles of incidence obtained from ROI-evaluation of MM-micromaps. Indices of fitted complex refractive indices  $N_\alpha$  and  $N_\beta$  refer to principal in-plane axes of index ellipsoid (c. figure left, inset in  $m_{34}$ -micrograph). Mind different y-scales for block diagonal and block off-diagonal MM-elements.

## Spectroscopic Imaging Mueller Matrix Ellipsometry


✓ acquire overall focused Mueller-Matrix micrographs  
 ✓ vary wavelength of illumination ( $\lambda=480-690$ nm in steps of 10nm)  
 ✓ micrographs reveal significant decrease of anisotropy for  $\lambda \rightarrow 680$  nm (c.  $m_{31}$  graphs below)  
 ✓ set ROI for stack-evaluation ("pixelshot")

ROI analysis


✓ ROI calculates average over the enclosed micrograph pixels for each wavelength and for each element of the measured normalized 3x4 Mueller Matrix  
 ✓ result: MM-SE spectrum (c. fig. below) of sample area of ~25  $\mu\text{m}^2$

ellipsometric modelling

✓ fitting in-plane complex refractive indices (wavelength by wavelength) yields optical dispersion in examined spectral range (c. fig. below). Note the convergence of  $\epsilon_2$  for  $\lambda > 600$ nm.  
 ✓ constant real value for out of plane refractive index  $N_z$  was used for modelling and calculation of in-plane values ( $N_z$  not accessible for single AOI).




Selected MM micrographs (elements  $m_{34}$  and  $m_{31}$ ) at different wavelengths of illumination. Green boxes highlight areas used for stack calculation of MM-spectrum.



Calculated in-plane dispersion functions, assuming constant  $N_z = 5.0 + i * 0$  and layer thickness of 43nm according to AFM data.


## AFM measurements

✓ atomic force microscopy (AFM) topography measurements yield layer thicknesses  $d$  of different parts of the flake  
 ✓ statistical analysis of central region (circular mask, c. image right):  
 $d = 43 \pm 3$  nm  
 ✓ line profiles yield approximate thicknesses of different plateaus



AFM measurements were carried out at University of Goettingen by Anna Sinterhauf and Dr. Martin Wenderoth.

Set mask for derivation of mean thickness and standard deviation of layer thickness in the flake's center.



Topography image obtained from AFM measurements (left) and horizontal line profiles with derived layer thicknesses (right).

## Summary

A microstructured flake of exfoliated black phosphorus was characterized by means of Imaging Mueller Matrix Ellipsometry (IMME). The flake-size was below the x-y-resolution limits of conventional ellipsometers. Thus, IMME is a technique to extend Mueller-matrix analysis of anisotropic 2D-materials to microstructured samples.

**IMME- $\theta$ -scans:**  
 ✓ non-zero off-diagonal MM-blocks of b-P-flake (MM-micromaps) reveal the anisotropy of thin-layer flake  
 ✓ recorded  $\theta$ -scans were successfully modeled by assuming a simple orthorhombic crystal structure with biaxial anisotropy for b-P-layer with one axis orthogonal to the sample's surface  
 ✓ in-plane complex refractive indices and orientations of in-plane optical axes were obtained for a single wavelength ( $\lambda=550$  nm) from fitting the model to IMME- $\theta$ -scans for three different angles of incidence using predetermined layer thickness (AFM) and real refractive index in z-direction  
 ✓ in-plane fit results are almost independent of refractive index in z-direction

**Spectroscopic IMME:**  
 ✓ once having found an ellipsometric model and in-plane orientations of optical axes obtained from IMME- $\theta$ -scans, a spectroscopic IMME measurement yielded the optical dispersions of in-plane components of b-P for the wavelength-range of 480-690 nm (assuming constant real refractive index for z-direction)  
 ✓ MM raw data and fitted dispersions both revealed significant reduction of in-plane anisotropy for spectral range of 600-700nm  
 ✓ Obtained dispersions might still depend on assumed value for refractive index in z-direction. Additional data at different AOI or sample rotations  $\theta$  might decouple in-plane and out-of-plane refractive indices upon numerical analysis.

References:  
 [1] A. Castellanos-Gomez, *J. Phys. Chem. Lett.*, vol.6, no.21 (Nov. 2015)  
 [2] A. Castellanos-Gomez et al., *2D Materials*, vol.1, no.2 (June 25, 2014)