


Brewster angle microscope
nanofilm_ultrabam



nanofilm_ultrabam with nanofilm_troughlarge

Brewster angle microscope

nanofilm_ultrabam


ABSTRACT

The nanofilm_ultrabam is the ultimate Brewster angle microscope designed for the air/liquid interface. It allows the direct visualization of Langmuir monolayers or adsorbate films. It also works on dielectric substrates like glass, quartz or similar materials.

The Ultrabam combines high resolution and overall focused real-time imaging. Advanced imaging optics provide fully focused images at 20-35 frames per second. A high performance camera and specific calibration algorithms enable quantitative measurements of reflectivity. Thus adsorption kinetics or thickness variations can be monitored.

The motorized analyzer offers the visualization of optical anisotropy, that is e.g. caused by long range molecular orientation order in monolayers.

The nanofilm_ultrabam is designed to match most troughs available. Regarding best performance and simple handling the nanofilm_trough/large is recommended. Vertical position changes of the instrument head are performed by a highly precise motorized lift.



nanofilm_ultrabam with trough, laser safety cabinet and support frame

APPLICATIONS


- Domains and order phenomena
- Transformation of monolayers into multilayered structures
- Photochemical reactions (e.g. photoisomerization)
- LB-films on solid structures
- Polymers and other materials that cannot be detected by fluorescence microscopy
- Adsorption kinetics
- In-situ polymerization in monolayers
- Phase separation
- Influences of subphase compositions (counter-ions) on monolayer structures
- Quality and homogeneity of LB-films
- ... and many more

FEATURES & BENEFITS


- Real-time overall focused images without focus scanner
- Direct visualization of the sample
- Lateral resolution down to 2 microns
- Imaging of anisotropic layers
- Designed for the air/liquid interface, but also works on dielectric substrates
- Variable angle-of-incidence ranging from 52° to 57°
- Large field of view
- On-line image processing
- AVI recorder
- Software interface for major trough systems
- Image files contain trough data
- Active vibration isolation system included




Graphical user interface of the nanofilm_ultrabam



Monolayer of 2-Monopalmitoleoyl-rac-glycerol at different surface pressure (left: $\pi = 3,06$ mN/m, middle: $\pi = 4,30$ mN/m, right: $\pi = 5,86$ mN/m)




Monolayer of 1-Monopalmitoleoyl-rac-glycerol at different surface pressure (left: $\pi = 2,68$ mN/m, middle: $\pi = 4,35$ mN/m, right: $\pi = 7,94$ mN/m)



QBAM tool used for quantitative Brewster angle microscopy studies




Calibration of the laser




Resolution of the imaging system according to USAF 1951 test target




Monolayer experiment with NIMA trough



Conventional imaging system used in standard Brewster angle microscopes



Sophisticated imaging system used in Ultra-bam



Geometrical correction

ACCESSORIES AND OPTIONS

- Laser safety cabinet
- Steel support frame
- nanofilm_trough
- Manual sample stage

Technical Specifications

nanofilm_ultrabam

CONFIGURATION

| | |
|---------------------------|--|
| TYPE: | Brewster angle microscope |
| ANGLE-OF-INCIDENCE RANGE: | 52° - 57° (equivalent to Brewster angle for substrates with $n = 1.3 \dots 1.54$), motorized, 0.001° resolution |
| LIGHT SOURCE: | 50 mW, 658 nm class IIIb broadband laser source for reduced interference fringes (other light sources on request) |
| IMAGING OPTICS: | True 2 micrometers resolution (horizontal image direction) according to Rayleigh's criterion; effective resolution of the camera 0.7 $\mu\text{m}/\text{pixel}$ (valid @ 53.1° AOI) |
| FIELD OF VIEW: | 720 x 400 μm |
| POLARIZER: | High quality Glan-Thompson prism, motorized, 0.001° resolution |
| ANALYZER: | Advanced high quality thin film polarizer, motorized, 0.001° resolution |
| CAMERA: | High performance CCD camera, 1360 x 1024 pixels, 2 x 2 binning mode for increased sensitivity, 20 fps (max: 35 fps @ 2 x 2 binning), 12 bits, GigE interface, variable exposure time and gain |
| MOUNTING: | Variable vertical position of the instrument head (alternative mounting options on request). Two-axis horizontal alignment of instrument head by fine-thread screws |
| Z-LIFT: | Precise motorized positioning of the instrument head, range 40 mm, adjustable lower end switch for safety reasons |
| VIBRATION ISOLATION: | halcyonics_variobasic_40 high performance active vibration isolation system (stable support structure and solid floor are included) |
| IMAGE PROCESSING: | Automatic background compensation, geometric image correction for unskewed images, contrast enhancement and image filtering, images contain trough data and BAM settings, AVI video recorder, including AccurionServer environment for simple access to measurement data |

ELECTRONICS

| | |
|-----------------|---|
| CONTROL: | nanofilm_ultrabam control software for simple access to all motorized components, camera etc.; Interface to KSVNIMA trough to import data on surface pressure etc.; Calibration routines for quantitative Brewster angle microscopy |
| ELECTRONICS | Based on nanofilm_ep4 electronics with embedded Linux controller |
| POWER SUPPLY: | 100 - 240 V, 50/60 Hz |
| PC AND MONITOR: | Up-to-date PC with Windows 7® operating system, pre-installed software - ready to use, LCD display 23" or larger |

EUROPE:

Accurion GmbH
Stresemannstrasse 30
37079 Goettingen, Germany

Phone: +49(0)551.999 60.0
Fax: +49(0)551.999 60.10
E-Mail: info@accurion.com
Web: www.accurion.com

NORTH AMERICA:

Accurion, Inc.
Inner Harbor Center
400 East Pratt Street, 8th floor
Baltimore MD 21202

Fax: +1-866.387.1714
E-Mail: info@accurion.com
Web: www.accurion.com

INDIA:

Accurion Scientific Instruments Pvt. Ltd.
Flat 402, WING 2, Balaji Lakeside Marvel
#1413/2A, /th Cross, 1st Main road
BHCS Layout, Uttarahalli
Bangalore 560 068, India

Phone: +91(0)98.450.04273
E-Mail: sharma@accurion.com
Web: www.accurion.com

CHINA:

Accurion Scientific Instruments (Shanghai) Co. Ltd.
Rm. 502, Xuhui Commercial Building,
No. 168, Yude Road
Xuhui District
Shanghai 200030, China

Phone: +86-(0)21.5017.9099
E-Mail: fz@accurion.cn
Web: www.accurion.com